



CPREzy™ CPR Assistance Device

#CPREZ01



#CPREZ01

The CPREzy CPR Assistance Device is unlike other manual External Chest Compression Assist Devices (ECCAD). Just turn it on, place it to the chest as suggested, and begin compressing. There is no waiting for voice prompts which divert your attention and waste precious rescue time. The CPREzy ensures you are giving correct compression depths regardless of the surface that your victim is laying upon.

Includes:

CPREzy Chest Pad & Battery

Black Belt Pouch

Instruction Leaflet

Compliant with current American Heart Association guidelines for CPR and external chest compressions.

FDA certified for training and rescue.



CPREzy™ Kit #CPREZ02

Includes:

Clam Shell case

CPREzy Chest Pad & Battery

CPREzy Mask & Mouthpiece

Gloves

Training CD/Instruction Leaflet

How it works

The Chest Pad has an audible chirp and visual light pacing system in the metronome component. The metronome speaker gives a clear chirp tone at 100 compressions-per-minute to match CPR administration guidelines. The audible chirp tone can be silenced while the visual light pacing system continues at the 100 compression-per-minute rate.

