



LIFE® SoftPac™

LIFE is Time™

**When minutes count
EMERGENCY OXYGEN
can be the difference between
LIFE OR DEATH.**

- ☒ **Light Weight...Only 6 Pounds!**
4 1/2" D x 7 1/2" W x 15" H
- ☒ Comfortable Handle & Shoulder Strap
- ☒ **Wall Mount Hook** for easy access
- ☒ **Constant reading supply gauge** through window (other units must be turned on to read gauge)
- ☒ **40 Minute Supply @ 6 Liters Per Minute**
- ☒ 250 Liter Aluminum "C" Cylinder
- ☒ **Cylinder is easily disposable/replaceable or low-cost refillable**
- ☒ Simple External Two-step instructions
Inhalator: administer supplemental oxygen for a breathing victim
Resuscitator: oxygen enriched mouth-to-mask CPR for a non-breathing victim
- ☒ **5 Year Warranty**
- ☒ **NON-PRESCRIPTION**



This unit could save your life.

Bag Color Choice:
White or Green

Make LIFE® SoftPac™ part of your first aid and safety program.



LIFE® CPR Mask
invertably fits both ADULT and CHILD (patented)

MODEL	FEATURES	Regulators
LIFE-2-612	6&12 LPM-"NORM"&"HIGH"(2 settings) provides to recommended 100% inspired O2	
LIFE-2-025	0 to 25 LPM Variable Flow (only for trained EMTs) 12 flow settings-0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25	
LIFE-2-101	Spare Cylinder Assembly complete with cylinder, valve, gauge & O2 fill	

LIFE® SoftPac™ MEETS OR EXCEEDS ALL GOVERNMENT AND FDA REQUIREMENTS: A NON-PRESCRIPTION DEVICE INTENDED FOR EMERGENCY USE.

"Knurled Knob" ON/OFF VALVE



Durable, Easy-to-Use, On-Off Valve
(With Constant Reading Gauge)

Uncompromisingly manufactured with quality materials, technology, and workmanship

FDA Regulation minimum capability is 6 LPM flow rate. LIFE® 612™ models deliver both the FDA minimum of 6 LPM. And the AHA recommended 100% inspired oxygen at 12 LPM. In just two simple settings which read "NORM" & "HIGH". "More is better (HIGH), but less lasts longer (NORM)". Only LIFE Corporation offers the simple 6&12 LPM flow rates.

American Heart Association Guidelines for CPR Oxygen

- ☒ "... give 100% inspired oxygen during Basic Life Support ... and Advanced Cardiac Life Support ... (BLS & ACLS) ... as soon as it becomes available."
- ☒ "... use supplemental oxygen...at a flow rate of 10-12 LPM ... to enable delivery of 100% inspired oxygen" (initially).