



Emergency Instruction Device  
2010 ECC/AHA Compliant!\*\*

Endorsed by the  
American  
Lifeguard  
Association



Approved by  
major training  
organizations.


# OSHA Best Practice: First Aid & CPR Training Drill Kit

#TRDK01

**To anyone responsible for training: stop struggling with scheduling, cost overruns, and retention risks.** First Voice has developed a system that will minimize your scheduling inconveniences, improve skill retention, AND reduce your training costs! Using our simple products, you can improve your responder programs, reduce your costs, and keep your responders at closer to 100% skills retention all year long.

Our Hands On CPR and First Aid Training system contains state-of-the-art products to ensure that your responders are always up-to-date and ready for action. Just 15-30 minutes every month guarantees that your responders know how to respond to any emergency.

Refreshing skills often with the Prestan Hands On CPR Feedback Manikin, First Voice Emergency Instruction Device (EID), the First Voice Training Scenarios Book, and an optional AED trainer upgrade in a convenient carrying case ensures that you are providing your responders with the best training and equipment at a combo, all-in-one price that won't break your budget.



Normal instructor CPR/First Aid class results.

- \* Increases skills retention
- \* Exceeds OSHA training requirements
- \* Proves "good faith" to OSHA inspectors

*The system comes with complimentary CPR First Aid or BBP skills assessment checklists for proper file documentation by your staff and management.*

*This system is a necessity to shore up OSHA hands-on training requirements when your organization does 100% online training and needs skills assessment verification and practice. It is a wonderful accessory to help your in-house responder team be well-prepared for any emergency.*



Comes with large duffel for all kit contents!



## Kit Includes:

### 1. The Emergency Instruction Device

The Emergency Instruction Device (EID #AVU5001x) is the heart of the First Voice products. At only 1.25 inches thick, the EID is backed by a 3-year warranty on all parts and is patent protected. All protocols are compliant with current, approved First Aid and ECC protocols, and this device is updateable. Utilizing yes/no questions quickly gets the trainee through practice training scenarios. It also provides: CPR metronome, AED refresher training, assessment reminders, basic First Aid refreshers, and other training drills.

### 2. The Training Scenarios Book

The Training Scenarios Book (#TS001) can be used as needed for a quick, 'anytime training' refresher, bringing the instructor to you. Covers medical, injury, and environmental emergencies and ensures skills practice for CPR and a majority of workplace first aid incidents. Comes with quizzes and answer key for your workplace documentation records.

### 3. The Prestan CPR Feedback Manikin

Realistic to the eye and the touch, the Prestan manikin (#PP-AM-100M) is unlike any other on the market. Prestan's revolutionary new CPR Rate Monitor allows for instant feedback to both instructor and student regarding the rate of chest compression and depth. This allows each student to gauge their rate of compressions on their own as well as allowing the instructor/trainer to monitor several students quickly and easily.

### 4. Skills checklist forms (CPR, First Aid, BBP)

### 5. Responder Handbook (V10570) 128 + Pages

### 6. CPR Kneeling Cushion (2 ea)

### 7. 10 Replacement lung bags for Prestan & Cleaning Instructions

### 8. 25 face shields, 100 cleaning wipes

### 9. Job assignment badges (for training drill roles.)

### 10. CPR/AED training drill example including Reviewer feedback scorecard

### 11. SAMPLE/DOTS Cue Cards

### 12. Responder CPR 2 pocket mask w/10pk replacement 1 way valves, 50pr nitrile gloves

\* According to the American Heart Association studies

\*\* American Heart Association 2010 Circulation Supplement Volume 122, Issue 18, Supplement 3, Nov. 2, 2010

- Part 5, S697: "... real-time CPR prompting and feedback technology such as visual and auditory prompting devices can improve the quality of CPR"

- Part 16, S920: "... nonresponders most frequently cited panic (37.5%) and fear of hurting the patient (9.1%) as reasons they were unable to perform."

- Part 16, S921: "even brief training increases the willingness to use an AED and improves performance."

- Part 16, S922: "Short video instruction combined with synchronous hands on practice is an effective alternative to instructor-led basic life support courses."

- Part 16, S923: "... ECC courses should be part of a larger continuing education and continuous quality improvement process."

- Part 16, S923: "Checklists or cognitive aids may be considered for use during actual resuscitation." " CPR prompt & feedback devices can be useful as part of an overall strategy to improve the quality of CPR during actual resuscitations."

## Recommended Add-on:

### Universal AED Trainer #GAEDT1

For refresher training on AEDs and sudden cardiac arrest (SCA) events at your organization. **Fits inside the Training System case!**

